


# Clean Swipe™ Test Coupons

## Description

For use in the Biotech and Pharmaceutical industry the Clean Swipe™ Test Coupon is critical when performing cleanability or swab recovery studies. The test coupon provides a representative sample of the material used for process equipment construction.

## Features

- Representative sample of material used on process equipment
- Serialized for tracking of coupon and identification
- 25 cm<sup>2</sup> Laser etched "Target Box" to indicate control area \*\*
- Individually bagged to keep surface protected and clean
- Raw Material Certificates provided with each order to document quality of material



Part Number	Material	Size	Polish Front	Polish Back	Target Box Laser Mark **
CPN-304L-25-32MP	Stainless Steel 304L	2.5" x 2.5" x 1/16" THK	32Ra MP	Unpolished	Yes
CPN-316L-25-32MP	Stainless Steel 316L	2.5" x 2.5" x 1/16" THK	32Ra MP	Unpolished	Yes
CPN-316L-25-20MP	Stainless Steel 316L	2.5" x 2.5" x 1/16" THK	20Ra EP	Unpolished	Yes
CPN-316L-25-20EP	Stainless Steel 316L	2.5" x 2.5" x 1/16" THK	20Ra EP	Unpolished	Yes
CPN-C22-25-32MP	Hastelloy C22	2.5" x 2.5" x 1/16" THK	32Ra MP	Unpolished	Yes
CPN-C276-25-32MP	Hastelloy C276	2.5" x 2.5" x 1/16" THK	32Ra MP	Unpolished	Yes
CPN-GLS-25-N	Borofloat™ Glass	2.5" x 2.5" x 3/16" THK	no polish	Unpolished	No
CPN-PTFE-25-N	Virgin PTFE	2.5" x 2.5" x 1/4" THK	no polish	Unpolished	No

\*\* - Target Box is a lasered outline of a 5cm x 5cm box that is marked on center of the coupon to aid with test evaluation

# TRU-FLOW

Your Process Partner

Broomfield, Colorado  
(877) 409-8590  
info@tru-flow.com  
www.tru-flow.com